

## 1. GENERAL INFORMATION

<b>IUPAC Name:</b>	2-(2,5-dimethoxy-4-propylphenyl)ethanamine
<b>CFR:</b>	1308.11; Schedule I
<b>CAS #:</b>	207740-22-5
<b>Synonyms:</b>	2,5-dimethoxy-4-propylphenethylamine
<b>Source:</b>	DEA Reference Material Collection
<b>Appearance:</b>	White powder (HCl)
<b>Kovat's Index:</b>	Pending
<b>UV<sub>max</sub>:</b>	289.9 nm

## 2. CHEMICAL AND PHYSICAL DATA

### 2.1 CHEMICAL DATA

Form	Chemical Formula	Molecular Weight	Melting Point (°C)
Base	C <sub>13</sub> H <sub>21</sub> NO <sub>2</sub>	223	Not Determined
HCl	C <sub>13</sub> H <sub>21</sub> NO <sub>2</sub> · HCl	259	213.3

### 3. ADDITIONAL RESOURCES

[Forendex](#)

[Wikipedia](#)

### 4. QUALITATIVE DATA

#### 4.1 NUCLEAR MAGNETIC RESONANCE


##### *Method NMR CDCl<sub>3</sub>*

*Sample Preparation:* Extract analyte with saturated sodium bicarbonate in D<sub>2</sub>O and deuteriochloroform (CDCl<sub>3</sub>) to ~10mg/mL.


**Instrument:** Varian Mercury 400 MHz NMR spectrometer with proton detection probe

**Parameters:** Spectral width: at least containing -3 ppm through 13 ppm  
Pulse angle: 90°  
Delay between pulses: 45 seconds  
Number of scans (NT): 8  
Number of steady state scans: 0  
Oversampling: 4 or more  
Shimming: automatic gradient shimming of Z1-4 shims  
Phasing, Drift Correction: automatic or manual


1H NMR: 2C-P HCl Lot # N1P32; CDCl<sub>3</sub>; 400MHz



1H NMR: 2C-P HCl Lot # N1P32; CDCl<sub>3</sub>; 400MHz



1H NMR: 2C-P HCl Lot # N1P32; CDCl<sub>3</sub>; 400MHz



## 4.2 GAS CHROMATOGRAPHY/MASS SPECTROMETRY

*Sample Preparation:* Dilute analyte to ~1 mg/mL base extracted in CHCl<sub>3</sub>.

**Instrument:** Agilent gas chromatograph operated in split mode with MS detector

**Column:** DB-1 MS; 30m x .25mm x .25µm

**Carrier Gas:** Helium at 1 mL/min

**Temperatures:** Injector: 280°C

MSD transfer line: 280°C

MS Source: 230°C

MS Quad: 150°C

Oven program:

1) 90°C initial temperature for 2.0 min

2) Ramp to 300°C at 14°C/min

3) Hold final temperature for 10.0 min

**Injection Parameters:** Split Ratio = 25:1, 1 µL injected

**MS Parameters:** Mass scan range: 34-550 amu


Threshold: 100

Tune file: stune.u

Acquisition mode: scan

**Retention Time:** 10.806 minutes


EI Mass Spectrum: 2C-P HCl Lot # N1P32



### 4.3 INFRARED SPECTROSCOPY (FTIR)

**Instrument:** FTIR with diamond ATR attachment (3 bounce)  
**Scan Parameters:** Number of scans: 32  
Number of background scans: 32  
Resolution:  $4\text{cm}^{-1}$   
Sample gain: 8  
Aperture: 150

FTIR ATR (Diamond, 3 Bounce): 2C-P HCl Lot # N1P32



FTIR ATR (Diamond, 3 Bounce): 2C-P HCl Lot # N1P32

