




## 5-Fluoro-MN-24

The Drug Enforcement Administration's Special Testing and Research Laboratory generated this monograph using structurally confirmed reference material.



### 1. GENERAL INFORMATION

**IUPAC Name:** 1-(5-fluoropentyl)-N-(naphthalen-1-yl)-1H-indole-3-carboxamide

**CAS#:** 1445580-60-8

**Synonyms:** MN-24F, CBM-2201, 5F-NNE1

**Source:** DEA Reference Material Collection

**Appearance:** White powder

**UV<sub>max</sub> (nm):** Not Determined

### 2. CHEMICAL AND PHYSICAL DATA

#### 2.1 CHEMICAL DATA

Form	Chemical Formula	Molecular Weight	Melting Point (°C)
Base	$C_{24}H_{23}FN_2O$	374	157.0



## 5-Fluoro-MN-24

The Drug Enforcement Administration's Special Testing and Research Laboratory generated this monograph using structurally confirmed reference material.




### 3. QUALITATIVE DATA

#### 3.1 NUCLEAR MAGNETIC RESONANCE

**Sample Preparation:** Dilute analyte to ~25 mg/mL in DMSO containing TMS for 0 ppm reference and maleic acid as quantitative internal standard.

**Instrument:** 400 MHz NMR spectrometer  
**Parameters:** Spectral width: at least containing -3 ppm through 13 ppm  
Pulse angle: 90°  
Delay between pulses: 45 seconds





## 5-Fluoro-MN-24

The Drug Enforcement Administration's Special Testing and Research Laboratory generated this monograph using structurally confirmed reference material.



### 3.2 GAS CHROMATOGRAPHY/MASS SPECTROMETRY


*Sample Preparation:* Dilute analyte ~4mg in 1mL Methanol

**Instrument:** Agilent gas chromatograph operated in split mode with MS detector  
**Column:** DB-1 MS (or equivalent); 30m x 0.25 mm x 0.25  $\mu$ m  
**Carrier Gas:** Helium at 1 mL/min  
**Temperatures:** Injector: 280°C  
MSD transfer line: 280°C  
MS Source: 230°C  
MS Quad: 150°C  
Oven program:  
1) 100°C initial temperature for 1.0 min  
2) Ramp to 300°C at 12 °C/min  
3) Hold final temperature for 25.0 min

**Injection Parameters:** Split Ratio = 20:1, 1  $\mu$ L injected  
**MS Parameters:** Mass scan range: 30-550 amu  
Threshold: 100  
Tune file: stune.u  
Acquisition mode: scan

**Retention Time:** 23.454 min

EI Mass Spectrum, 5F-MN-24 Lot # RM-131115-02





## 5-Fluoro-MN-24

The Drug Enforcement Administration's Special Testing and Research Laboratory generated this monograph using structurally confirmed reference material.





### 3.3 INFRARED SPECTROSCOPY (FTIR)

**Instrument:** FTIR with diamond ATR attachment (3 bounce)

**Scan Parameters:**  
Number of scans: 32  
Number of background scans: 32  
Resolution: 4 cm<sup>-1</sup>  
Sample gain: 8  
Aperture: 150

FTIR, ATR (Diamond, 3 Bounce), 5F-MN-24 Lot# RM-131115-02





## 5-Fluoro-MN-24

The Drug Enforcement Administration's Special Testing and Research Laboratory generated this monograph using structurally confirmed reference material.



### 4. ADDITIONAL RESOURCES

Shevyrin, V.; Melkozerov, V.; Nevero, A.; Eltsoz, O.; Shafran, Y.; Analytical characterization of some synthetic cannabinoids, derivatives of indole-3-carboxylic acid. **Forensic Science International**. 2013, p1-10