

## 1. GENERAL INFORMATION

<b>IUPAC Name:</b>	4-methylhexan-2-amine
<b>CFR:</b>	Not Scheduled (4/2013)
<b>CAS #:</b>	105-41-9
<b>Synonyms:</b>	1,3-dimethylamylamine, methylhexanamine, 2-amino-4-methylhexane, 1,3-dimethylpentylamine, 4-methyl-2-hexylamine
<b>Source:</b>	DEA Reference Material Collection
<b>Appearance:</b>	White powder (HCl)
<b>Kovat's Index:</b>	Pending
<b>UV<sub>max</sub>:</b>	Not Determined

## 2. CHEMICAL AND PHYSICAL DATA

### 2.1 CHEMICAL DATA

Form	Chemical Formula	Molecular Weight	Melting Point (°C)
Base	C <sub>7</sub> H <sub>17</sub> N	115	Not Determined
HCl	C <sub>7</sub> H <sub>17</sub> N · HCl	151	122.9

### 3. ADDITIONAL RESOURCES

[Wikipedia](#)

### 4. QUALITATIVE DATA

#### 4.1 NUCLEAR MAGNETIC RESONANCE


##### *Method NMR Dimethylfumarate/DMSO*

*Sample Preparation:* Dilute analyte to ~10 mg/mL in DMSO containing TMS for 0 ppm reference and dimethylfumarate as quantitative internal standard.


**Instrument:** 400 MHz NMR spectrometer

**Parameters:** Spectral width: at least containing -3 ppm through 13 ppm  
Pulse angle: 90°  
Delay between pulses: 45 seconds

1H NMR: DMAA HCl Lot N1P7; DMSO; 400MHz





1H NMR: DMAA HCl Lot N1P7; DMSO; 400MHz



***NMR Analytical Observation***

DMAA has two chiral carbons; therefore, two diastereomers are possible. Diastereomers produce slightly different chemical shifts for proton and carbon. Both diastereomers are present in the above spectra.



## 4.2 GAS CHROMATOGRAPHY/MASS SPECTROMETRY

*Sample Preparation:* Dilute analyte to ~1 mg/mL base extracted into chloroform

***Instrument:*** Agilent gas chromatograph operated in split mode with MS detector

***Column:*** DB-1 MS or equivalent; 30m x .25mm x .25 $\mu$ m

***Carrier Gas:*** Helium at 1 mL/min

***Temperatures:*** Injector: 280°C

MSD transfer line: 280°C

MS Source: 230°C

MS Quad: 150°C

Oven program:

1) 100°C initial temperature for 1.0 min

2) Ramp to 300°C at 12°C/min

3) Hold final temperature for 9.0 min

***Injection Parameters:*** Split Ratio = 20:1, 1  $\mu$ L injected

***MS Parameters:*** Mass scan range: 30-550 amu


Threshold: 100

Tune file: stune.u

Acquisition mode: scan

***Retention Time:*** 2.100 minutes


EI Mass Spectrum: DMAA HCl Lot N1P7



### 4.3 INFRARED SPECTROSCOPY (FTIR)

**Instrument:** FTIR with diamond ATR attachment (3 bounce)  
**Scan Parameters:** Number of scans: 32  
Number of background scans: 32  
Resolution:  $4\text{cm}^{-1}$   
Sample gain: 8  
Aperture: 150

FTIR ATR (Diamond, 3 Bounce): DMAA HCl Lot N1P7



FTIR ATR (Diamond, 3 Bounce): DMAA HCl Lot N1P7

